

LNLS script and Py4Syn



What is the LNLS script?

- An easy way to start / stop IOCs
- An easy way to access IOCs EPICS shell
- An easy way to know the status of all IOCs running in the machine
- Based on procServ
 - No need to know the mapping between IOCs and telnet ports
 - LNLS script takes care of port numbers transparently

- Start VMware
- Start the VM for the tutorial
- If VMware asks, tell it you did a copy of the VM
- User = guest / password = guest
- If needed, user = root / password = root

What is Py4Syn?

Py4Syn (*Python For Synchrotron*) is:

- A Python package
- Used as a high level tier
- A Device manipulation tool
- Used for scans
- Used to see charts in real-time
- Simple

Motivation

Requirements for ideal scenario:

- Free.
- Open source.
- Easy to maintain and improve.
- Standardized but scalable to attend the most different needs.
- Support to beamline's devices.
- Motion control, counters, scans, pseudo-motors.
- Realtime charts.
- Simple syntax, clear and most closer to what the users are used to as possible.

- **Devices**
 - Motors, Counters, CCDs, etc...
- **Utility functions**
 - Motor movement
 - Counting
 - Scans
 - Charts
 - Statistics
- **Simple, direct and user focused documentation.**



```
from py4syn.utils.motor import *

createMotor('mtop', 'SOL:DMC1:m3') # Create a motor named mtop

umv('mtop', 10) # Move mtop to absolute position 10
                # Wait for arrival
                # Show position while moving

wa() # Show the position of all motors created

currentPosition = wmr('mtop') # Read the current position

mv('mtop', 20, wait=False) # Move mtop to absolute position 20
                           # Do not wait for arrival
                           # Do not show position while moving

# Here code can run while motor is moving
```

Complete list in:

http://py4syn.readthedocs.org/en/latest/utils/utils_motor.html

Pseudo-Motors

Virtual motors that are described by equations relating one or more motors.

Note: There is no automatic velocity and acceleration settings to make all motors arrive in the target at the same time. The total movement time is the time given by the motor that arrives last.

Pseudo-Motor structure:

- Name
- Description
- Readback equation (how the pseudo position is calculated)
- Target equation dictionary (how to define the target of each axis that is part of the pseudo-motor)

Example: create a pseudo-motor to control gap and offset of a slit

```
from py4syn.utils.motor import *
```

```
createMotor('mtop', 'SOL:DMC1:m3') # Create motor mtop
```

```
createMotor('mbot', 'SOL:DMC1:m4') # Create motor mbot
```

```
gapRBV = 'A[mtop]-A[mbot]' # Equation to describe gap size
```

```
gapTopTarget = 'A[offset]+T[gap]/2.0' # Equation to define mtop target
```

```
gapBotTarget = 'A[offset]-T[gap]/2.0' # Equation to define mbot target
```

```
offsetRBV = 'A[mtop]-A[gap]/2.0' # Equation to describe offset position
```

```
offsetTopTarget = 'T[offset]+A[gap]/2.0' # Equation to define mtop target
```

```
offsetBotTarget = 'T[offset]-A[gap]/2.0' # Equation to define mbot target
```


```
createPseudoMotor('offset', 'off slit 1', offsetRBV, {mtop: offsetTopTarget, mbot: offsetBotTarget})
```

```
createPseudoMotor('gap', 'gap slit 1', gapRBV, {mtop: gapTopTarget, mbot: gapBotTarget})
```

```
wa() # Show motor position for all motors, including pseudo
```

```
umv('gap', 5) # Move gap to 5 and wait
```

```
wa() # Show motor position for all motors, including pseudo
```



Counters

```
from py4syn.epics.ScalerClass import *
from py4syn.utils.counter import *

scaler = Scaler('SOL:SCALER', 10, 'scaler1') # Create a scaler

createCounter('det', scaler, 3) # Create a counter named det
createCounter('mon', scaler, 10) # Create a counter named mon

ct(1)# Run a 1 preset time count in all counters
      # Show results at the end

counting = ctr(1000, use_monitor=True)# Run a counting and return data

disableCounter('det') # Disable det counter
```

Complete list: http://py4syn.readthedocs.org/en/latest/utils/utils_counter.html

There are 3 scanning types available:

- **Scan**

- In a Scan, the number of points is the same for all used devices.
- Devices will search for its targets at the same time.
- There are no restrictions to the number of devices or points used.
- There are no restrictions to devices types being used as long as the IScannable interface is implemented.

- **Mesh**

- In a mesh the number of points is independent for each used device.
- For each point of a device, the next device will do a complete scan before the former device goes to the next point.
- There are no restrictions to the number of devices or points used.
- There are no restrictions to devices types being used as long as the IScannable interface is implemented.

- **Timescan**

- In a timescan no devices but counters are used. This scan runs in time defined intervals.
- It is possible to define a delay time between measurements.

scan function

Standard format of scan function:

```
scan(dev1, start1, end1, dev2, start2, end2,... ,devn, startn, endn, points, time)
```

Where:

devN is any IScannable device.

startN is the initial value.

endN is the final value.

points is the number of intervals between points to be done (Py4Syn includes initial and final values).

time is the integration time to be configured in used counters.

Optionally, the pair **startN**, **endN** can be changed by an array, as long as all devices use the same number of points and each array size equals the parameter “points”.

The same is applied to **time**, where an array of integration time can be used, as long as the array size equals the parameter “points”.

In a scan, the total number of points is given by “points” + 1.

scan function

Example:

- `scan('m1', 0, 2, 'm2', 2, 6, 2, 0.1)`

Where:

- m1 will go from 0 to 2 with steps of value 1 $((2-0)/2)$
- m2 will go from 2 to 6 with steps of value 2 $((6-2)/2)$
- 3 points in total (2 points + 1 final point)
- 0.1 seconds of integration time

Results in:

Point	m1	m2
0	0.0	2.0
1	1.0	4.0
2	2.0	6.0

Standard format of mesh function:

`mesh(dev1, start1, end1, points1, dev2, start2, end2, points2,... ,devN, startN, endN, pointsN, time)`

Where:

devN is any IScannable device.

startN is the initial value.

endN is the final value.

pointsN is the number of intervals between points to be done (Py4Syn includes initial and final values).

time is the integration time to be configured in used counters.

Optionally, the pair **startN, endN** can be changed by an array and, in this case, the parameter **pointsN** is omitted.

In mesh, the total number of points is given by $(points1+1)*(points2+1)*(...)*(pointsN+1)$.

mesh function

Example:

- `mesh('m1', 0, 2, 2, 'm2', 2, 6, 2, 0.1)`

Where:

- m1 will go from 0 to 2 with steps of value 1 $((2-0)/2)$
- m2 will go from 2 to 6 with steps of value 2 $((6-2)/2)$
- 9 points in total (2 points + 1 final point) to m1 and m2, and for mesh, total = $p_{m1} * p_{m2} = 3 * 3$
- 0.1 seconds of integration time

Results in:

Point	m1	m2
0	0.0	2.0
1	0.0	4.0
2	0.0	6.0
3	1.0	2.0
4	1.0	4.0
5	1.0	6.0
6	2.0	2.0
7	2.0	4.0
8	2.0	6.0

Standard format of timescan function:

`timescan(time, delay, repeat)`


Where:

time is the integration time to be configured in used counters. Default value is 1 second.

delay is the time to wait until the next counting. Default value is 1 second.

repeat is the number of counting to be done. Default value is infinite.

Devices classification



Complete list: <http://py4syn.readthedocs.org/en/latest/epics.html>

In Py4Syn library, devices are divided into two groups:

IScannable

The IScannable interface defines some standard operations that all devices intended to be used during a scan must implement. For example:

- **setValue** - Define the target value
- **getValue** - Read the current value
- **wait** - Wait until the target is hit
- **getLowLimit** - Read low limit
- **getHighLimit** - Read high limit

ICountable

The ICountable interface defines some standard operations that all devices intended to be used as counters must implement:

- **getValue** – Read the current value
- **setCountTime** - Define integration time
- **setPresetValue** - Define preset value
- **startCount** – Start a counting
- **stopCount**- Stop a counting
- **canMonitor** – Define if it can be used as a monitor
- **canStopCount** - Define if counting can be interrupted
- **isCounting** – Verify if counting is being done
- **wait** – Wait until counting ends

These special interfaces allow a generalization of devices for scanning functions independently of its type: a motor, a furnace, a MCA, a CCD, etc.

Independently of the scan type, many customizations are possible to the actions by means of user defined functions (*Callbacks*).

There are 7 available callbacks for all scan types:

- **Pre Scan Callback**
 - Called before the start of a scan.
- **Pre Point Callback**
 - Called before any device is commanded to its target.
- **Pre Operation Callback**
 - Called after all devices arrive in the target and before the start of the counters.
- **Operation Callback**
 - Called while the counters are counting.
- **Post Operation Callback**
 - Called after the counting and before the chart updates.
- **Post Point Callback**
 - Called after all operations and before the start of a new point.
- **Post Scan Callback**
 - Called after the end of the scan.

Example: Defining a callback to be called during the counting

```
from py4syn.utils.scan import *

def myCallback(**kwargs): # Function that will be called during the operation.
    scanObject = kwargs['scan'] # Reference to scan object.
    indexArray = kwargs['idx'] # Reference to the point index.
    positionArray = kwargs['pos'] # Reference to the devices position.
    print('Callback message')

setPreScanCallback(defaultPreScanCallback) # Configuring PreScanCallback to default.

setOperationCallback(mycallback) # Configuring OperationCallback to our function.

scan('m1', 0, 180, 10, 1) # Running the scan with the configuration.
```

There is a need to show data in real-time and the need that data rendering would not generate lateness in running scans.

- Real-time rendering.
- Many charts in the same windows.
- Data overlap in the same chart.

Example 1: Simple chart


```
import time
from py4syn.utils.plotter import *

pl = Plotter("Teste de Plotter", daemon=False) # Create a Plotter (Window)
pl.createAxis("Plot 1", xlabel="Eixo X", ylabel="Eixo Y", grid=True,
line_style="-", line_marker="o", line_color="blue", label="Dados") # Create a
new axis

for i in range(0, 20):
    pl.plot(i, i*i, axis=1) # Put a new point in chart
    time.sleep(0.1)
```

Complete list: http://py4syn.readthedocs.org/en/latest/utils/utils_plotter.html

Charts



Example 2: Advanced chart

```
import time
from py4syn.utils.plotter import *

N = 100 # Number of points in chart
pl = Plotter('Plot Sample e I0', daemon=True) # Create a Plotter for Sample and I0
pl.createAxis(title="Energy Scan", xlabel="Energy", ylabel="I0", grid=True, line_style="--",
line_marker="x", line_color="blue", label="I0") # Create an axis for I0
pl.createAxis(title="", xlabel="Energy", ylabel="Sample", grid=True, line_style=":", line_marker="o",
line_color="black", label="Sample") # Create an axis for Sample



pl1 = Plotter('Plot Norm', daemon=True) # Create a Plotter for normalized data
pl1.createAxis(title="", xlabel="Energy", ylabel="Norm", grid=True, label="1")# Create an axis for
Norm

idx = 1
for i in range(0, 10): # Number of scans we will run
    if(i >= 1):
        idx += 1
        pl1.createAxis(title="", parent=1, label=str(i+1)) # For each scan, create a new axis for Norm

pl.clear(1) # For each scan, clean I0 chart
pl.clear(2) # For each scan, clean Sample chart

for ii in range(N): # For each point
    i0 = random.random()*100 # Generate a value for I0
    sample = random.random()*127 # Generate a value for Sample
    norm = (sample/i0)/10 # Calculate normalization
    pl.plot(ii, i0, axis=1) # Update I0 chart
    pl.plot(ii, sample, axis=2) # Update Sample chart
    pl1.plot(ii, norm, axis=idx) # Update Norm chart in the correct axis
    time.sleep(0.01)
```

Charts



Special Variables

In scan module it is possible to access special parameters by its read and write methods. For example:

- SCAN_DATA – Dictionary that saves all scan data.
- SCAN_CMD – Last scan command run.
- XFIELD / YFIELD – X/Y axis of the automatic chart.
- PEAK / PEAK_AT – Value and position of the peaks among the data.
- FIT_SCAN – Define if data might or might not be automatically fit at the scan's end.
- PLOT_GRAPH – Define if data might or might not be automatically shown in chart.
- among others...

Complete list: http://py4syn.readthedocs.org/en/latest/utils/utils_scan.html#special-variables-in-the-scan-module

Focusing in better performance and the lowest possible dead-time, data is kept in memory during data registration. Only at the end of a scanning (by error or success) data is transferred to disk. To change this default behavior, one must use the function `setPartialWrite(True)`. Doing so, data is saved at each iteration.

The standard file format is the same used by SPEC and PyMCA because beamlines in LNLS are historically used to it. This format is supported by many analysis programs, too.

Example:

```
#E 1413896783
#D Tue Oct 21 11:06:23 2014
#C py4syn User = user.name
#C This is a comment, ignore it
#S 1 scan(gap, 10, 11, 10, 0.01)
#D Tue Oct 21 11:06:15 2014
#T 0.01 (Seconds)
#N 7
#L points gap seconds I0 cyber testeField timestamp
0 10.0 7.495 0.50 0.733 N/A 2014-10-21 11:06:17.423802
1 10.1 5.973 0.40 0.496 N/A 2014-10-21 11:06:17.988520
...
...
```

Saving data using a different format:

- Create a class that implements FileWriter interface
- Methods to be implemented:
 - `writeHeader(self)` – Write the header to a file
 - `writeData(self, partial=False, idx=-1)` – Write data to file
 - `close(self)` – Close the file
- Pass the new class to scan object:
`scan.setFileWriter(writer)`

`writeData` must use data information provided by `SCAN_DATA` dictionary. Its structure is described in http://py4syn.readthedocs.org/en/latest/utils/utils_scan.html#handling-the-data

Currently the interpolation module has only a Gaussian one and gives information such:

- Full width at half maximum / Position
- Maximum / Position
- Minimum / Position
- Center of mass / Position

Example:

```
import pylab as pl
from py4syn.utils.fit import *
data = 'testdata.txt'
X = pl.loadtxt(data); # Load a data file
x = X[:,0]; # Read X values
y = X[:,7]; # Read Y values

pkv, pkp, minv, minp, fwhm, fwhmp, com = fitGauss(x, y) # Run fitGauss
print("Peak ", pkv, " at ", pkp) # Print values and its calculated positions
print("Min ", minv, " at ", minp)
print("Fwhm ", fwhm, " at ", fwhmp)
print("COM = ", com)
```

Py4Syn

Py4Syn is open-source. Code can be get at
<https://github.com/hhslepicka/py4syn>

Complete documentation: <http://py4syn.lnls.br>

Published paper:

<http://scripts.iucr.org/cgi-bin/paper?S1600577515013715>

